Disposable Inspiratory One-way Mouthpiece with

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






7

e.g. achieved = 55

Target
Clinically
Effective
Flow Rate
30-90 L/Min

8

Caution

Consult 'Instructions for Use'.

This product complies with the essential requirements of the medical devices directive (93/42/EEC).

en

Instructions for Use

Introduction

The In-Check DIAL G16 kit is an inhalation airflow training meter system that can help educate and assess patients who use inhaler devices.

Inhaler devices are designed to deliver medication to the respiratory tract, and the speed of inhalation through them (the inspiratory flow) can have a significant effect on the quantity of drug delivered and the clinical efficacy of the product.

The In-Check DIAL G16 kit is designed to simulate the "internal resistance" of common inhaler devices, and assess inspiratory flow. These assessments enable the healthcare professional to encourage patients to modify their inspiratory technique (by inhaling with more, or less effort), in order to achieve a flow rate consistent with clinical efficacy. The coloured 'flow' icons show the clinically effective flow ranges for each different inhaler device. These 'flow' icons do not imply any comparison between devices.

Patients that cannot achieve the suggested inspiratory flow for their inhaler may not gain maximum benefit from their prescribed medication, and healthcare professionals may wish to take this factor into account when selecting with the patient, the device that is the most suitable.

Inspiratory Flow and Clinically Effective Flow Range

The inspiratory flow through an inhaler is one of the factors that will influence the clinical effect of the drug delivery from that device. The most effective delivery occurs when the patient achieves a flow within the clinically effective flow range. Flow rates outside this range, may result in a diminished deposition and clinical efficacy.

For further information, see: www.haag-streit.com/clement-clarke

In-Check DIAL G16

is an assessment device that has a selectable resistance from high to low, shown by the coloured 'flow' icons on the back cover of this booklet, calibrated to enable the assessment of airflow as if the patient was using certain different inhalers.

The "conventional pMDI"  position can be used to approximate the resistance when a conventional pMDI is attached to a holding chamber or spacer device.

When used with the enclosed mouthpiece, In-Check DIAL G16 provides an assessment of inspiratory flow through five classes of device resistance. However, the resistance is not precise for all inhalers, e.g. the resistance of the Genuair is slightly lower than the resistance of the Turbohaler Mk III, but they are both medium resistance inhalers. To mimic the precise resistance of a particular inhaler, specific adaptors are available to use with the DIAL set to 'conventional pMDI'.

Should new types of inhaler become available, then In-Check DIAL G16 kit can be used to assess the new inspiratory flow, by knowing the resistance of the new inhaler.

IMPORTANT

As with any inhalation device, it is important to check for loose foreign objects before the device is used to prevent debris entering the mouth/throat. The transparent material used in the construction of the In-Check DIAL G16 enables the user to make a visual check before inhalation. Patients should be prevented from exhaling through the device prior to use to avoid damaging the device.

To Reset the In-Check DIAL G16

Hold the instrument vertically with the mouthpiece uppermost, so that the rounded end of the meter can be tapped firmly against the other hand or a horizontal surface, such as a table.

This will dislodge the magnetic resetting weight, which will return the red cursor to the start position. The device must now be turned through 180 degrees to return the magnetic weight to its resting position.

Do Not attempt to reset the In-Check DIAL G16 by shaking/tapping at any other orientation.

- this action will cause serious damage to the piston and pointer and end the life of the instrument.

How to use the In-Check DIAL G16 kit - see diagrams

1. Reset the In-Check DIAL G16. See '**To Reset the In-Check DIAL G16**' section.
2. Align the dial selector with the desired coloured icon – an audible "click" should be heard.
3. Attach a clean mouthpiece. Disposable one-way inspiratory mouthpieces are preferred.
4. Ask the patient to exhale fully.
5. Ask the patient to seal their lips around the mouthpiece. According to the inhaler chosen, instruct the patient to inhale in the manner recommended by the inhaler manufacturer.
6. Record the inspiratory flow from the position of the red cursor against the scale. Reset, and repeat two more times, ensuring correct technique each time.
7. Compare values achieved with target flows for that device. To operate an inhaler device correctly, the patient should be able to achieve a flow rate within the clinically effective range. *e.g. Accuhaler: Clinically Effective Flow Rate 30-90 L/Min.
8. If after repeated training the patient is not able to achieve these values, then the healthcare professional may wish to assess the patient's ability to use an alternative type of inhaler.

Performance Accuracy

Accuracy +/- 10% or 10 L/min (whichever is the greater) and repeatability of +/- 5 L/min.

Cleaning your In-Check DIAL G16

As the In-Check device is designed for inspiratory flow assessment, the risk of cross-infection is greatly reduced when compared to an expiratory peak flow meter because clean air is being sucked in. However, the meter can become compromised if a patient accidentally blows or coughs into the device.

In an effort to reduce the risk of cross-contamination to its lowest possible level we recommend the use of the inspiratory one-way valve mouthpieces or Spirometry type anti-bacterial/anti-viral filters.

Cleaning Methodology

For personal use your In-Check device can be cleaned by immersing the device in warm (but not hot) mild detergent solution for 2-3 minutes (maximum 5 minutes). Agitate the meter to ensure thorough cleaning. Rinse in warm water and shake to remove any excess water. It is important to rinse thoroughly to prevent salt spots appearing on the inside of the body and the spindle.

To shake excess water from the In-Check device, hold only at the end furthest away from the device selector.

Allow to dry thoroughly before using again.

For more intensive cleaning/disinfection between patients the In-Check device can also be cleaned and disinfected using "Steris ProKlenz NpH" or "Rely+On/PeraSafe".

The In-Check device can be sterilized using the "Low Temperature hydrogen peroxide plasma - Sterilization System".

It is important to follow the manufacturer's instructions and recommended contact times for the above products.


Cleaning Frequency

The following recommendations for disinfection frequencies are presented as a guide only. In practice, the person responsible for the clinical wellbeing of the patient should consider the specific circumstances of the next patient and the risk of infection posed by cross-contamination.

If used with a one-way valve mouthpiece or anti-bacterial filter the In-Check device should be cleaned weekly.

If used with an open mouthpiece the device should be cleaned between patients.

Interpreting information from the In-Check DIAL G16

In-Check DIAL G16 was modified in January 2016, in response to customer feedback, to provide information about a wide range of new inhaler devices and the inspiratory flow rates that are associated with clinical efficacy. The new scale is shown above.

The information provided is for guidance only. It does not imply that any particular product will be clinically effective. There are other aspects than inspiratory flow that contribute towards clinical efficacy. The information is based on published clinical studies. In some cases studies may have been performed on a particular formulation in a specific device and may be suggestive that other formulations will behave similarly.

DPIs

Laboratory studies demonstrate that DPIs have a degree of flow-dependence, meaning that delivery may be affected by increasing flow rates. This does not always equate to improved clinical efficacy. A pale arrow represents the intention that flow rates of a higher level may be beneficial for DPIs.


pMDIs

The pMDIs have been represented in a single block from 20-60 L/Min because the clinical studies support this. The pale arrow is pointed in the opposite direction for pMDIs, indicating that a slower rather than higher flow rate is beneficial in this type of inhaler. The devices are individually represented on the back cover.

Adaptors (specific device mimics)

These adaptors listed below can be used with the In-Check DIAL G16 :

Turbohaler Mk II	- T	- 3109781
HandiHaler	- H	- 3109859
Diskhaler 4 Blister	- D4	- 3109788
Clickhaler	- C	- 3109786
Genuair	- G	- 3109861
Breezhaler	- B/F	- 3109783
Easyhaler (M) Monotherapy	- EH	- 3109787
Accuhaler	- A/A/D	- 3109822

Set DIAL to the pMDI to use with device specific adaptor.

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Capstick, T. G., & Clifton, I. J. (2014). Inhaler technique and training in people with chronic obstructive pulmonary disease and asthma. *Expert review of respiratory medicine*.

Banno M, Ibata H, Niimi T, Sato S, Matsushita R. The Usefulness of Inspiratory Flow Rate during Inhalation Corticosteroid Therapy in Asthma. *Respiration* 2009;78:387-392.

Janssens W, vanden Brande P, Hardeman E, et al. Inspiratory flow rates at different levels of resistance in elderly COPD patients. *Eur Respir J*. 2008;31:78-83.

Zhang QL, Zheng JP, Pan WH, He H, Chen GL, An JY, Yuan LF, Luo DF. Peak inspiratory flow generated through different analogue dry powder inhalers in Shenzhen healthy preschool children. *Zhonghua Er Ke Za Zhi*. 2008 Feb;46(2): 98-103.

Raid A.M. R.A.M. Al-Showair, Walid Y. W.Y. Tarsin, Khaled H. K.H. Assi, Stanley B. S.B. Pearson, Henry Chrystyn. Can all patients with COPD use the correct inhalation flow with all inhalers and does training help? *Respiratory Medicine* - November 2007 Vol. 101, Issue 11, Pages 2395-2401.

Walid Y. Tarsin, Stanley B. Pearson, Khaled H. Assi, Henry Chrystyn. Emitted dose estimates from Seretide(R) Diskus(R) and Symbicort(R) Turbuhaler(R) following inhalation by severe asthmatics. *International Journal of Pharmaceutics*, Volume 316, Issues 1-2, 19 June 2006, Pages 131-137.

K.N. Desager, J. Geldhof, R. Claes, W. De Backer. Measurement of Inspiratory Flow Through Three Different Dry Powder Inhalation Devices Using In-Check™ in Children with Asthma. *Pediatric Asthma, Allergy & Immunology*. March 2006, Vol. 19, No. 1: 6-13.

Andrea S. Melani; Letizia S. Bracci; Marcello Rossi. Reduced Peak Inspiratory Effort through the Diskus and the Turbuhaler due to Mishandling is Common in Clinical Practice. *Clin Drug Invest*. 2005;25(8):543-549.

Amirav I, Newhouse MT, Mansour Y. Measurement of peak inspiratory flow with In-Check DIAL device to simulate low-resistance (Diskus) and high-resistance (Turbuhaler) dry powder inhalers in children with asthma. *Pediatr Pulmonol*. 2005 May;39(5):447-51.

Y. Adachi, Y.S. Adachi, T. Itazawa, J. Yamamoto, G. Murakami, T. Miyawaki.

Measurement of peak inspiratory flow rates with an In-Check™ meter to identify preschool children's ability to use dry powder inhalers; Diskus and Turbuhaler. *The Journal of Allergy and Clinical Immunology* - February 2004 (Vol. 113, Issue 2, Supplement, Page S114.

Please Note

All trademarks and product names are the property of their respective owners.

Handihaler is a registered trademark of **Boehringer Ingelheim Pharma GmbH & Co. KG**.

Easyhaler is a registered trademark of **Orion Corporation**.

NEXThaler is a registered trademark of **Chiesi Farmaceutici S.P.A.**

Twisthaler is a registered trademark of **MSD International Holdings GmbH**.

Turbuhaler and **Turbohaler** are registered trademarks of **AstraZeneca AB**.

Clickhaler is a trademark of **Tianjin Kinnovata Pharmaceutical Company Limited**.

Spiromax and **Respiclick** are registered trademarks of **IVAX International B.V.**

Genuair, **Novolizer** and **Pressair** are registered trademarks of **AstraZeneca AB**.

Ellipta is a registered trademark of **GlaxoSmithKline group of companies**.

Diskhaler is a registered trademark of **GlaxoSmithKline group of companies**.

Accuhaler and **Discus** are registered trademarks of **GlaxoSmithKline group of companies**.

Breezhaler and **Aerolizer** are registered trademarks of **Novartis AG**.

Repimat is a registered trademark of **Boehringer Ingelheim Pharma GmbH & Co. KG**.

Alliance Tech Medical
PO Box 6024 / Granbury, TX 76049
Tele: 1.800.848.8923 / 817.326.6357
Fax: 817.326.2182
www.alliancetechmedical.com

©Copyright 2019 Clement Clarke International Ltd.

Part no. 3109305 Issue 9 Sep 2019

Icon	Product	Icon	Product	Device Resistance
	*Handihaler®		Forspiro®	
	*Easyhaler®		Ellipta®	
	NEXThaler®		TurboSpin®	
	Twisthaler®		*Diskhaler®(4)	
	*Turbohaler® P		*Accuhaler®	
	Turbohaler® S		*Breezhaler®	
	*Clickhaler™		k-haler®	
	Spiromax®		Respirat®	
	*GenuAir®			
Please Note: All trademarks and product names are the property of their respective owners, see IFU booklet for details. *Indicates device specific adaptor is available, see IFU booklet.				

Assessing inspiratory flow rate for clinical efficacy:

Select appropriate resistance setting, inhale through meter, assess achieved flow rate.


For DPIs values between 30-90 L/min are generally associated with clinical efficacy.

For pMDIs values between 20-60 L/min are preferred.


For Handihaler 20 L/min is associated with clinical efficacy.

Bibliography:

1. Pedersen S, Dubus JC, Crompton G. 5) Inhaler selection in children with asthma. *Prim Care Respir J* 2010;**19** (3):209-16
2. Lavorini F, Levy ML, Corrigan C, Crompton G. 6) Training tools for inhalation devices. *Prim Care Respir J* 2010;**19** (4):335-41
3. Capstick, T. G., & Clifton, I. J. (2014). Inhaler technique and training in people with chronic obstructive pulmonary disease and asthma. *Expert review of respiratory medicine*.
4. Sanders, M. J. (2017). Guiding Inspiratory Flow: Development of the In-Check DIAL G16, a Tool for Improving Inhaler Technique. *Pulmonary medicine*, 2017.
5. Mahler, D. A. (2017). Peak inspiratory flow rate as a criterion for dry powder inhaler use in chronic obstructive pulmonary disease. *Annals of the American Thoracic Society*, 14(7), 1103-1107.
6. Laube, B. L., Janssens, H.M., de Jongh, F.H., Devadason, S.G., Dhand, R., Diot, P., ... & Chrystyn, H. (2011). What the pulmonary specialist should know about the new inhalation therapies. *European Respiratory Journal* **37** (6) 1308-1417.

Inhaler Resistance Range

Alliance Tech Medical
PO Box 6024, Granbury, TX 76049
1.800.848.8923 / 817.326.6357 Fax: 817.326.2182
www.alliancetekmedical.com

Alliance
TECH MEDICAL