

In-Check™

The In-Check is a portable inspiratory flow meter that monitors both the disease and the response to treatment. It is designed for use in:

- Measurement of oral inspiratory flow rate (peak and sub-maximal) in order to assess suitability for inhaled medication
- Measurement of nasal inspiratory flow for assessing nasal airway patency
- Assessment of response to nasal provocation
- Home monitoring of Peak Nasal Inspiratory Flow (PNIF)



Disposable Mouthpieces

Economical and convenient. Reduce risk of cross-infection, helping to protect patient and equipment of contamination. One-way valve mouthpieces are available for inspiratory use.



Filters

Provide protection against possible cross-infection helping to protect patient and equipment.

Masks

Can be added to measuring devices

Ordering Information

	Order #	Quantity	NHRIC#
In-Check	3109750	Oral - each	3 8462109750 3
	3109751	Nasal with Mask	3 8462109751 0
Disposable Mouthpieces	3122200	Adult - 100/bag	3 8462122200 4
	3123200	Pediatric - 100/bag	3 8462123200 3
	3122069	One-way insp. - 200/bx	3 8462122069 7
Filters	5551100	Filter - 100/case	10 38462551100 2
Masks	2005100	Adult Mask - each	3 8462005100 1
Disposable Adapters	2004100	Pediatric Mask - each	3 8462004100 2
	3122120	Pediatric/Small Mouth Adult - bag of 20	3 8462122120 5
	3122220	Pediatric/Small Mouth Adult - box of 100	3 8462122220 2



PO Box 6024 / Granbury, TX 76049
 1.800.848.8923 / 817.326.6357
 Fax: 817.326.2182
 www.alliancetechnical.com

To place all orders
 CALL: 1.800.848.8923 / 817.326.6357

Accuracy

Each In-Check is individually calibrated to ensure a high degree of accuracy (+/- 10% or 10 L/min - whichever is greater). Inspiratory flow rates in the range of 30 L/min - 370 L/min can be measured.

US FDA 510k approved. Conforms to international standards for peak flow meter performance.

Features

- High-quality medical-grade stainless steel and plastic construction
- Sterilizable plastic mouthpiece accepts different mask sizes
- Simple enough for home patient use, rugged enough for routine clinic use.

Cleaning

The In-Check should be cleaned by washing in warm soapy water, rinsed and dried thoroughly.

Call today for Alliance Tech Medical's complete line of COPD, respiratory, asthma management products, sanitizers, splash shields, ventlab chambers and many other items.

Distributed By: