

In-Check™ Flo-Tone

The pMDI is a widely used popular inhaler format; yet most patients do not use their pMDIs correctly. Many inhale too rapidly, some cannot co-ordinate inhalation and actuation easily.

The need for patients to use inhalers correctly is becoming widely recognised.

Product Information

Designed to improve patient use of pMDIs. In-Check Flo-Tone is a simple, innovative training device that makes it easy for healthcare professionals to help patients to achieve correct flow rate and master inhale/actuate co-ordination.

With In-Check Flo-Tone attached to the pMDI, the patient's gentle inhalation causes a whistle to sound, this is a signal to actuate the inhaler. As the patient maintains the whistle the aerosol is drawn deep into the lungs.

What is this Product used for?

The Flo-Tone device is designed to be used by healthcare professionals to coach patients, helping them improve their inhaler technique; by signalling the correct flow rate for actuation and guiding them to sustain their inhalation.



Product Features and Benefits

- For use with placebo pMDI training inhalers
- Patient learns to inhale at correct flow rate
- Inhaler training can be conducted quickly
- Whistle provides clear positive signal to actuate
- Low-cost, simple yet effective method of advanced inhaler technique training
- Currently fits many pMDIs, other sizes to follow

Product Suitability



Technical Specifications

Measurement Range: 0.02 - 0.24 kPa

Weight: 5g

Accuracy: 100% tested

Part Code

In-Check Flo-Tone (pack of 10): 3152002

Peak Inspiratory Flow