



In-Check™



## The first multi-purpose, mechanical inspiratory flow measurement device for clinical and in-home use.

**Easy to use, clean and reset, the In-Check allows accurate measurement of both oral and nasal inspiratory flow rates that help to assess respiratory condition and treatment.**



**Improve your technique**

- For use by children and adults
- Allows home monitoring of Peak Nasal Inspiratory Flow (PNIF)
- Transparent body permits visual inspection
- Accepts standard or disposable mouthpieces along with one-way valved disposable mouthpieces
- Optional mask available for measuring nasal inspiratory flow

Help patients maintain their respiratory health with the In-Check inspiratory flow meter.



# In-Check™

The In-Check is a portable inspiratory flow meter that monitors both the disease and the response to treatment. It is designed for use in:

- Measurement of oral inspiratory flow rate (peak and sub-maximal) in order to assess suitability for inhaled medication
- Measurement of nasal inspiratory flow for assessing nasal airway patency
- Assessment of response to nasal provocation
- Home monitoring of Peak Nasal Inspiratory Flow (PNIF)



## Disposable Mouthpieces

Economical and convenient. Reduce risk of cross-infection, helping to protect patient and equipment of contamination. One-way valve mouthpieces are available for inspiratory use.



## Filters

Provide protection against possible cross-infection helping to protect patient and equipment.

## Masks

Can be added to measuring devices

## Ordering Information

	Order #	Quantity	NHRIC#
<b>In-Check</b>	3109750	Oral - each	3 8462109750 3
	3109751	Nasal with Mask	3 8462109751 0
<b>Disposable Mouthpieces</b>	3122200	Adult - 100/bag	3 8462122200 4
	3123200	Pediatric - 100/bag	3 8462123200 3
	3122069	One-way insp. - 200/bx	3 8462122069 7
<b>Filters</b>	5551100	Filter - 100/case	10 38462551100 2
<b>Masks</b>	2005100	Adult Mask - each	3 8462005100 1
<b>Disposable Adapters</b>	2004100	Pediatric Mask - each	3 8462004100 2
	3122120	Pediatric/Small Mouth Adult - bag of 20	3 8462122120 5
	3122220	Pediatric/Small Mouth Adult - box of 100	3 8462122220 2

## Accuracy

Each In-Check is individually calibrated to ensure a high degree of accuracy (+/- 10% or 10 L/min - whichever is greater). Inspiratory flow rates in the range of 30 L/min - 370 L/min can be measured.

US FDA 510k approved. Conforms to international standards for peak flow meter performance.

## Features

- High-quality medical-grade stainless steel and plastic construction
- Sterilizable plastic mouthpiece accepts different mask sizes
- Simple enough for home patient use, rugged enough for routine clinic use.

## Cleaning

The In-Check should be cleaned by washing in warm soapy water, rinsed and dried thoroughly.

**Call today for Alliance Tech Medical's complete line of COPD, respiratory, asthma management products, sanitizers, splash shields, ventlab chambers and many other items.**



**Alliance Tech Medical, Inc.™**

PO Box 6024 / Granbury, TX 76049  
 1.800.848.8923 / 817.326.6357  
 Fax: 817.326.2182 / www.alliancetechmedical.com

**Eastern Office** 302.838.2638 **Central Office** 630.910.0993  
**Southern Office** 817.326.3183 **Western Office** 800.848.8923

Distributed By: