



Easy on-PC

This easy-to-operate spirometer uses the power of your PC, laptop, or tablet, leveraging premium ultrasound technology for a complete spirometry solution.



Features

Added TrueFlow™ Advantages

- Not influenced by heat, humidity, barometric pressure, or altitude
- No cleaning, no maintenance
- Extremely high accuracy for low flows

Essential Features

- Flexible EMR capabilities
- Inspiratory and expiratory real-time curve(s)
- Challenge testing with manitol and methacholine protocols
- Pediatric incentives
- Trending
- Selectable predicted values and interpretation
- Automated quality control
- Quick test
- Multilingual



Incentive Screens

System Requirements

Operating system:	Windows XP SP3, Windows Vista, Windows 7, Windows 8 (32 and 64 bit)*
Processor:	Intel-/Pentium or equivalent, 1GHz
Interface:	USB 2.0 or earlier
Screen resolution:	min 1024 x 600
Hard disk capacity:	Installation/System 1GB, Memory 4GB
RAM:	min 512 MB - 1GB recommended

* Only "full" Windows 8 version is supported by our software for tablets.

Systems

Part Numbers	Description
ATM2700-3 Easy on-PC spirometer	Spirometry sensor, software, carry case
ATM2700-3K Easy on-PC spirometer kit	Spirometry sensor, software, carry case, (2) cases of 50 spirettes, (2) bags of 25 nose clips, 3-Liter Calibration Syringe, 3 year extended warranty

Disposables

Part Numbers	Description
ATM2050-1 spirettes	Case of 50
ATM2050-5 spirettes	Case of 200

Test Modes

- FVC, FVL, SVL, MVV, post testing
- Tidal breathing selectable

Proven Worldwide

- NASA, COPDGene, BOLD
- Meets ATS/ERS, SSA Disability and NIOSH/OSHA standards.



Learn more about Easy on-PC at alliancetechnical.com

PO Box 6024, Granbury, TX 76049 | 1.800.848.8923 / 817.326.6357 | Fax: 817.326.2182