

# Pulmonary Function Filters



AllFLOW™ Filters provide variety, ease-of-use and affordability. Available with round or **oval** patient side mouth ends



Introducing  
AllFLOW™ Oval

AllFLOW™ filters meet or exceed the highest standards in PFT filtering with minimal resistance and low dead space.

Our Port Color Variety (PCV)™ reduces contamination risk with filter mounting. The distinguishable color provides correct connection with equipment the first time. Quality spirometry and lung function tests are accurately obtained with 8 individual port sizes.

The AllFLOW™ filter medium traps bacteria and viruses efficiently with electrostatic properties.

- Low Resistance
- Minimal Dead Space
- Unique filter design accommodates Small Mouth Adapter for adults and children without compromising filter function
- 99.99% bacterial and viral filtration
- 8 standard port sizes defined by color designed to fit a wide selection of equipment.
- Filters exceed ATS guidelines for all criteria

Connect with us online




[www.alliancetekmedical.com](http://www.alliancetekmedical.com)

# Ordering Information

Port Side Color	Manufacturer Cross Filter	OD	ID	Quantity	Shape	Order#	NHRIC#/NDC
Grey	InCheck™, InCheck DIAL™, Multispiro™, CreativeBiomedics, Vitalograph, Schiller, MicroDirect, SDI (Pulmogard II™), Most CosMed Models, MIR (All Models), P.K. Morgan 232, Future Med (Spirobank), PB-400, PB-600, CDX Spiro 850, Spirometrics (Flowmate, Flowmate+, SMI I, Vacumed (Micro, Micro Plus, Micro GP-GL, Micro Loop, Pony)	30	26.3	100/box	Round Oval	5551100 7771100	10 38462551100 2 10 03846277110 4
Red	Brentwood, Burdick, and most Fukuda Models, Spirometrics (some models)	43.5	40	100/box	Round	5552100	10 38462552100 1
Orange	Koko, KoKo Legend, Keystone, Koko-Trek, CyberMedic Full Gas Systems, Collins 4G and 3G Full Gas Systems.	48.5	45	100/box	Round Oval	5553100 7773100	10 38462553100 0 10 38462773000 8
Medium Blue	SensorMedics, 2200,2400, 2450, 2800, others, Including Vmax, Keystone I & II, Pail, Cybermedics Moose, Gould 5004	34	30	100/box	Round Oval	5554100 7774100	10 38462554100 9 10 038462774105
Medium Green	Jaeger Models	33.5	30	100/box	Round Oval	5555100 7775100	10 38462555100 8 10 03846277510 2
Yellow	Collins DC-2	28.3*	28.7	100/box	Round Oval	5556100 7776100	10 38462556100 7 10 03846277610 9
Bright Green	Jaeger Master Screen Series, New	32.5	29.7	100/box	Round Oval	5558100 7778100	10 38462558100 5 10 03846277810 3
Blue-Grey	MedGraphics preVent™	37.95	34.60 (with 34.8 undercut)	100/box	Round Oval	5559100 7779100	10 38462559100 6 10 038462779100

\* 2 Step Port Side, Stepped up to 34.8 for Collins Equipment



Port Side Color	Filter Kit	Quantity	Order#	NHRIC#/NDC
Grey	1000 filter with noseclip & mouthpiece	50/box	5551595	50 38462551595 4
Orange	3000 filter with noseclip & mouthpiece	50/box	5553595	50 38462553595 2
Medium Blue	4000 filter with noseclip & mouthpiece	50/box	5554595	50 38462554595 1
Medium Green	5000 filter with noseclip & mouthpiece	50/box	5555595	50 38462555595 0
Yellow	6000 filter with noseclip & mouthpiece	50/box	5556595	50 38462556595 9
Bright Green	8000 filter with noseclip & mouthpiece	50/box	5558595	50 38462558595 7

**All Filter Kits Contain Filter, Nose Clip & Disposable Latex Free Mouthpiece**

## Specifications:

### Bacterial & Viral


**Filtration Efficiency:** 99.99%

**Dead Space:** Approx. 35ml (cc)

**Resistance:** 0.23 cm H<sub>2</sub>O/L/sec approx.

**Material:** K resin case  
Electrostat® filtrate medium

Description	Quantity	Order#	NHRIC#/NDC
Padded Nose Clip	100/bag	5550047	3 8462550047 4
Latex Free Rubber Disposable Mouthpiece	50/box	5550050	3 8462550050 1
Small Mouth Adapter	100/box	3122220	3 8462122220 2
Filter-Fit Disposable Cardboard Mouthpiece	100/bag	3125500	



Distributed By:



Alliance Tech Medical, Inc.

PO Box 6024 / Granbury, TX 76049

1.800.848.8923 / 817.326.6357

Fax: 817.326.2182 / www.alliancetekmedical.com

Eastern Office 302.838.2638 Central Office 800.848.8923

Southern Office 817.326.3183 Western Office 800.848.8923