

# Pulmonary Function Filters

## ALLFLOW™



**Alliance**  
TECH MEDICAL



ALLFLOW Filters provide variety, ease-of-use and affordability. Available with round or oval patient side mouth ends.

**ALLFLOW filters meet or exceed the highest standards in PFT filtering with minimal resistance and low dead space.**



Our Port Color Variety (PCV)™ reduces contamination risk with filter mounting. The distinguishable color provides correct connection with equipment the first time. Quality spirometry and lung function tests are accurately obtained with 8 individual port sizes. The ALLFLOW filter medium traps bacteria and viruses efficiently with electrostatic properties.


- Low Resistance
- Minimal Dead Space
- Unique filter design accommodates Small Mouth Adapter for adults and children without compromising filter function
- 99.99% bacterial and viral filtration
- 8 standard port sizes defined by color designed to fit a wide selection of equipment
- Filters exceed ATS guidelines for all criteria

# Ordering Information

## Filters

Port Side Color	Manufacturer Cross Filter	OD	ID	Quantity	Shape	Order#
Grey	InCheck™, InCheck DIAL™, Multispiro™, CreativeBiomedics, Vitalograph, Schiller, MicroDirect, SDI (Pulmogard II™), Most CosMed Models, MIR (All Models), P.K. Morgan 232, Future Med (Spirobank), PB-400, PB-600, CDX Spiro 850, Spirometrics (Flowmate, Flowmate+, SMI I, Vacumed) (Micro, Micro Plus, Micro GP-GL, Micro Loop, Pony)	30mm	26.3mm	100/box	Round Oval	5551100 7771100
Orange	Koko, KoKo Legend, Keystone, Koko-Trek, CyberMedic Full Gas Systems, Collins 4G and 3G Full Gas Systems	48mm	45mm	100/box	Round Oval	5553100 7773100
Medium Blue	SensorMedics, 2200,2400, 2450, 2800, others, Including Vmax, Keystone I & II, Pail, Cybermedics Moose, Gould 5004	34mm	30mm	100/box	Round Oval	5554100 7774100
Medium Green	Jaeger Models	33.5mm	30mm	100/box	Round Oval	5555100 7775100
Yellow	Collins DC-2	28.3*mm	28.7mm	100/box	Round Oval	5556100 7776100
Dark Blue	Morgan/Vitalograph Models		28.6mm	100/box	Round Oval	5557100 7777100
Bright Green	Jaeger Master Screen Series, New	32.5mm	29.7mm	100/box	Round Oval	5558100 7778100
Blue-Grey	MedGraphics preVent™	37.95mm	34.60mm (with 34.8mm undercut)	100/box	Round Oval	5559100 7779100

\* 2 Step Port Side, Stepped up to 34.8 for Collins Equipment



## All Filter Kits Contain Filter, Nose Clip & Disposable Latex Free Mouthpiece


Port Side Color	Filter Kit	Quantity	Order #
Grey	1100 filter with noseclip & mouthpiece	50/box	5551595
Orange	3100 filter with noseclip & mouthpiece	50/box	5553595
Medium Blue	4100 filter with noseclip & mouthpiece	50/box	5554595
Medium Green	5100 filter with noseclip & mouthpiece	50/box	5555595
Yellow	6100 filter with noseclip & mouthpiece	50/box	5556595
Dark Blue	7110 filter with noseclip & mouthpiece	50/box	5557595
Bright Green	8100 filter with noseclip & mouthpiece	50/box	5558595
Blue-Grey	9100 filter with noseclip & mouthpiece	50/box	5559595

## Accessories:

Description	Quantity	Order#
Padded Nose Clip	100/bag	5550047
Latex Free Rubber Disposable Mouthpiece	50/box	5550050
Small Mouth Adapter	100/box	3122220
Filter-Fit Disposable Cardboard Mouthpiece	100/bag	3125500

## Specifications:

Bacterial & Viral Filtration Efficiency:	99.99%
Dead Space:	Approx. 35ml (cc)
Resistance:	0.23 cm H <sub>2</sub> O/L/sec approx.
Material:	Technostat filtrate medium



Distributed by:



PO Box 6024 / Granbury, TX 76049  
 1.800.848.8923 / 817.326.6357  
 Fax: 817.326.2182  
 www.alliancetechnical.com