

# Revital-Ox™ RESERT™ R60 Solutions Test Strips

## All it takes is a minute: literally.



The Revital-Ox® RESERT® High Level Disinfectant has a **NEW AND IMPROVED** test solution strip.

### Easier to Use

Use on solutions at 20° C/68° F or higher. No need to cool solution before testing

### Cost Effective

12 month shelf life provides inventory flexibility

### Compliant

Dip before you soak. Test Revital Ox RESERT High Level Disinfectant solutions before each use.

Stay compliant to the Association for the Advancement of Medical Instrumentation and the Society of Gastroenterology Nurses and US disinfectant guidelines.

### FAST!

In exactly one minute, the test strip will tell you if the solution passes or fails.



**PASS**  
**FAIL**



### Order Now

Product Code	Description	Quantity
PCC053	Revital-Ox RESERT R60 Solution Test Strips	1 bottle containing 60 test strips
PCC054	Revital-Ox RESERT R60 Solution Test Strips	2 bottles containing 60 test strips/bottle

For more information, contact your **STERIS** representative or visit [steris.com/resertr60](http://steris.com/resertr60).

The Revital-Ox™ RESERT™ R60 Test Strips are part of an unbroken chain of innovative products and solutions from STERIS.

That is the power of...

